


EDGEONE™ FIRE™



PERFORMANCE, FLEXIBILITY & STRENGTH

EdgeOne Fire™ features our Heat-Treated Fire-Wire™ NiTi and is designed to shape canals in a reverse-reciprocating motion. EdgeOne Fire works with existing handpieces operated at the WaveOne® Gold parameter. EdgeOne Fire™ is available in sizes Small, Primary, Medium, and Large. Files have a parallelogram shaped cross section with a variable taper. There are three files in each pack with the option of an assortment pack or single sized pack. Available in lengths 21mm, 25mm and 31mm.

PERFORMANCE, FLEXIBILITY & STRENGTH



EdgeOne Fire™ - Sterile 3 Pack

Size	21mm	25mm	31mm	Gutta Percha (Pack of 60)	Paper Points (Pack of 60)
20 Small	EN-SEOFSM21	EN-SEOFSM25	EN-SEOFSM31	EN-EOFSMGP	EN-EOFSMPP
25 Primary	EN-SEOFPRI21	EN-SEOFPRI25	EN-SEOFPRI31	EN-EOFPRI GP	EN-EOFPRI PP
35 Medium	EN-SEOFMED21	EN-SEOFMED25	EN-SEOFMED31	EN-EOFMEDGP	EN-EOFMEDPP
45 Large	EN-SEOFLG21	EN-SEOFLG25	EN-SEOFLG31	EN-EOFLG GP	EN-EOFLG PP

EdgeOne Fire™ GlidePath - Non Sterile 4 Pack

Size	21mm	25mm	31mm
19	EN-EOGP1921	EN-EOGP1925	EN-EOGP1931

EdgeOne Fire™ Assorted - Sterile 3 Pack

Size	21mm	25mm	31mm
SM, PRI, MED	EN-SEOFAST21	EN-SEOFAST25	EN-SEOFAST31

