

EDGESEQUEL SAPPHIRE™



PERFORMANCE, FLEXIBILITY & STRENGTH

EdgeSequel Sapphire™ has a triangular cross section with an open flute design and is a constant taper up to the cutting flutes (cutting flute length is 15 mm). Comes in .04 15-90 and .06 15-50 taper sizes and tips. Available in lengths 21mm, 25mm and 31mm with four same size files in each pack. EdgeSequel Sapphire will help elevate your technique and practice.

PERFORMANCE, FLEXIBILITY & STRENGTH

EdgeSequel Sapphire™ - Sterile 4 Pack

ISO Size / Taper	21mm	25mm	31mm	Gutta Percha (Pack of 60)	Paper Points (Pack of 60)
15/.04	EN-SES150421	EN-SES150425	EN-SES150431	EN-ESGP1504	EN-ESPP1504
20/.04	EN-SES200421	EN-SES200425	EN-SES200431	EN-ESGP2004	EN-ESPP2004
25/.04	EN-SES250421	EN-SES250425	EN-SES250431	EN-ESGP2504	EN-ESPP2504
30/.04	EN-SES300421	EN-SES300425	EN-SES300431	EN-ESGP3004	EN-ESPP3004
35/.04	EN-SES350421	EN-SES350425	EN-SES350431	EN-ESGP3504	EN-ESPP3504
40/.04	EN-SES400421	EN-SES400425	EN-SES400431	EN-ESGP4004	EN-ESPP4004
45/.04	EN-SES450421	EN-SES450425	EN-SES450431	EN-ESGP4504	EN-ESPP4504
50/.04	EN-SES500421	EN-SES500425	EN-SES500431	EN-ESGP5004	EN-ESPP5004
55/.04	EN-SES550421	EN-SES550425	EN-SES550431	EN-ESGP5504	EN-ESPP5504
60/.04	EN-SES600421	EN-SES600425	EN-SES600431	EN-ESGP6004	EN-ESPP6004
70/.04	EN-SES700421	EN-SES700425	EN-SES700431	EN-ESGP7004	EN-ESPP7004
80/.04	EN-SES800421	EN-SES800425	EN-SES800431	EN-ESGP8004	EN-ESPP8004
90/.04	EN-SES900421	EN-SES900425	EN-SES900431		
15/.06	EN-SES150621	EN-SES150625	EN-SES150631	EN-ESGP1506	EN-ESPP1506
20/.06	EN-SES200621	EN-SES200625	EN-SES200631	EN-ESGP2006	EN-ESPP2006
25/.06	EN-SES250621	EN-SES250625	EN-SES250631	EN-ESGP2506	EN-ESPP2506
30/.06	EN-SES300621	EN-SES300625	EN-SES300631	EN-ESGP3006	EN-ESPP3006
35/.06	EN-SES350621	EN-SES350625	EN-SES350631	EN-ESGP3506	EN-ESPP3506
40/.06	EN-SES400621	EN-SES400625	EN-SES400631	EN-ESGP4006	EN-ESPP4006
45/.06	EN-SES450621	EN-SES450625	EN-SES450631	EN-ESGP4506	EN-ESPP4506
50/.06	EN-SES500621	EN-SES500625	EN-SES500631	EN-ESGP5006	EN-ESPP5006

EdgeSequel Sapphire™ Assorted - Sterile 4 Pack

Size	21mm	25mm	31mm
Assorted .04/15, 20, 25, 30"	EN-SESSAST0421	EN-SESSAST0425	EN-SESSAST0431
Assorted .04/35, 40, 45, 50"	EN-SESSAST0421L	EN-SESSAST0425L	EN-SESSAST0431L
Assorted .06/15, 20, 25, 30"	EN-SESSAST0621	EN-SESSAST0625	EN-SESSAST0631
Assorted .06/35, 40, 45, 50"	EN-SESSAST0621L	EN-SESSAST0625L	EN-SESSAST0631L

