


EDGETAPER PLATINUM™



PERFORMANCE, FLEXIBILITY & STRENGTH







EdgeTaper Platinum™ features our heat-treated Firewire NiTi, which gives amazing flexibility, capable of navigating 90 degree curves and has no bounce back to preserve apical anatomy. It has incredible strength and resistance to cyclic fatigue.

EdgeTaper Platinum™ file system has a bloated triangular cross section with a progressive changing taper. Available in lengths 21mm, 25mm, 31mm with six of the same size file in each pack. Experience the new standard in file safety and canal centering ability.

PERFORMANCE, FLEXIBILITY & STRENGTH



EdgeTaper Platinum™ - Sterile 6 Pack

Size	ISO Size	19mm	21mm	25mm	31mm	Gutta Percha (Pack of 60)	Paper Points (Pack of 60)
SX	19	EN-SETSX19HT					
S1 	17		EN-SETS121HT	EN-SETS125HT	EN-SETS131HT		
S2	20		EN-SETS221HT	EN-SETS225HT	EN-SETS231HT		
F1 	20		EN-SETF121HT	EN-SETF125HT	EN-SETF131HT	EN-EEGPF1	EN-EEPPF1
F2 	25		EN-SETF221HT	EN-SETF225HT	EN-SETF231HT	EN-EEGPF2	EN-EEPPF2
F3 	30		EN-SETF321HT	EN-SETF325HT	EN-SETF331HT	EN-EEGPF3	EN-EEPPF3
F4 	40		EN-SETF421HT	EN-SETF425HT	EN-SETF431HT	EN-EEGPF4	EN-EEPPF4
F5 	50		EN-SETF521HT	EN-SETF525HT	EN-SETF531HT	EN-EEGPF5	EN-EEPPF5
Assorted SX, S1, S2, F1, F2, F3			EN-SETP21MM	EN-SETP25MM	EN-SETP31MM		

