


GUTTA PERCHA POINTS - Pack of 60

Family	Tip Size	Taper Size	SKU
EdgeTaper™ & EdgeTaper Platinum™			
F1	20		EN-EEGPF1
F2	25		EN-EEGPF2
F3	30		EN-EEGPF3
F4	40		EN-EEGPF4
F5	50		EN-EEGPF5
EdgeFile™ .04			
	15	04	EN-SGP1504
	20	04	EN-SGP2004
	25	04	EN-SGP2504
	30	04	EN-SGP3004
	35	04	EN-SGP3504
	40	04	EN-SGP4004
	45	04	EN-SGP4504
	50	04	EN-SGP5004
EdgeFile™ .06			
	15	06	EN-SGP1506
	20	06	EN-SGP2006
	25	06	EN-SGP2506
	30	06	EN-SGP3006
	35	06	EN-SGP3506
	40	06	EN-SGP4006
	45	06	EN-SGP4506
	50	06	EN-SGP5006
EdgeOne Fire™			
	Small		EN-EOFSMGP
	Primary		EN-EOFPRIGP
	Medium		EN-EOFMEDGP
	Large		EN-EOFLGPP
EdgeSequel Sapphire™			
20	04		EN-ESGP2004
	06		EN-ESGP2006
25	04		EN-ESGP2504
	06		EN-ESGP2506
30	04		EN-ESGP3004
	06		EN-ESGP3006
35	04		EN-ESGP3504
	06		EN-ESGP3506
40	04		EN-ESGP4004
	06		EN-ESGP4006
45	04		EN-ESGP4504
	06		EN-ESGP4506
50	04		EN-ESGP5004
	06		EN-ESGP5006

PAPER POINTS - Pack of 60

Family	Tip Size	Taper Size	SKU
EdgeTaper™ & EdgeTaper Platinum™			
F1	20		EN-EEPPF1
F2	25		EN-EEPPF2
F3	30		EN-EEPPF3
F4	40		EN-EEPPF4
F5	50		EN-EEPPF5
EdgeFile™ .04			
	15	04	EN-SPP1504
	20	04	EN-SPP2004
	25	04	EN-SPP2504
	30	04	EN-SPP3004
	35	04	EN-SPP3504
	40	04	EN-SPP4004
	45	04	EN-SPP4504
	50	04	EN-SPP5004
EdgeFile™ .06			
	15	06	EN-SPP1506
	20	06	EN-SPP2006
	25	06	EN-SPP2506
	30	06	EN-SPP3006
	35	06	EN-SPP3506
	40	06	EN-SPP4006
	45	06	EN-SPP4506
	50	06	EN-SPP5006
EdgeOne Fire™			
	Small		EN-EOFSMPP
	Primary		EN-EOFPRIPP
	Medium		EN-EOFMEDPP
	Large		EN-EOFLGPP
EdgeSequel Sapphire™			
20	04		EN-ESPP2004
	06		EN-ESPP2006
25	04		EN-ESPP2504
	06		EN-ESPP2506
30	04		EN-ESPP3004
	06		EN-ESPP3006
35	04		EN-ESPP3504
	06		EN-ESPP3506
40	04		EN-ESPP4004
	06		EN-ESPP4006
45	04		EN-ESPP4504
	06		EN-ESPP4506
50	04		EN-ESPP5004
	06		EN-ESPP5006

THE NEW GOLD STANDARD FOR PERFORMANCE, FLEXIBILITY & STRENGTH



EdgeFile X7™

STERILE - PACK OF 6

- Heat-treated FireWire™ NiTi gives amazing flexibility, capable of 90° curves
- 2-8 times more cyclic fatigue resistant than other file systems.*
- Tips 17 - 45 .04 & .06 tapers
- Available in 21mm, 25mm & 29mm
- Max 1mm flute diameter
- Unmatched flexibility
- No shape memory
- Universal compatibility



EdgeGlidePath™

STERILE - PACK OF 4

- Heat-treated FireWire™ NiTi
- Replacing as many as 8 instruments to create a faster glidepath
- Replaces #15, 20, 25, 30, 35 hand files in a step-back technique and a #1, 2 & 3 Gates Glidden
- Available in 21mm, 25mm, and 31mm



EdgeOne Fire™

STERILE - PACK OF 3

- Heat-treated FireWire™ NiTi
- EdgeOne Fire™ is designed to shape canals in a reciprocating motion, in a similar way to leading reciprocating files brands
- EdgeOne Fire™ can be used with your existing motor and hand piece reciprocating settings
- Excellent flexibility, capable of 90° curves
- No bounce back to preserve canal anatomy
- Sizes: Small (20), Primary (25), Medium (35) & Large (45)
- Available in 21mm, 25mm and 31mm



EdgeTaper Platinum™

STERILE - PACK OF 6

- Heat-treated FireWire™ NiTi gives amazing flexibility, capable of 90° curves
- Familiar technique
- 2X more resistance to cyclic fatigue compared to competitor files*
- SX, S1, S2, F1, F2, F3, F4, F5
- 19mm, 21mm, 25mm, 31mm



EdgeSequel Sapphire™

STERILE - PACK OF 4

- Heat-treated FireWire™ NiTi
- 2X more resistant to cyclic fatigue compared to competitor files*
- Triangular cross section with open flute design
- Constant Taper
- ISO Sizes 15-90 .04 tapers
- ISO Sizes 15-50 .06 tapers
- Available in 21mm, 25mm and 31mm



EdgeV-Taper HT™

NON STERILE - PACK OF 6

- Heat-treated FireWire™ NiTi
- The Original Minimal preparation file
- Changing Variable Taper
- Tips 20, 22, 25, 30, 35, 40, 45 and 50
- GlidePath 13/.03, 14/.03 and 17/.04
- 21mm, 25mm and 31mm (29mm)

