

GUTTA PERCHA POINTS - Pack of 60

| Family | Tip Size | Taper Size | SKU |
|---|----------|------------|---------------|
| EdgeTaper™ & EdgeTaper Platinum™ | | | |
| F1 | 20 | | EN-ETGPF1-1 |
| F2 | 25 | | EN-ETGPF2-1 |
| F3 | 30 | | EN-ETGPF3-1 |
| F4 | 40 | | EN-ETGPF4-1 |
| F5 | 50 | | EN-ETGPF5-1 |
| EdgeFile™ .04 | | | |
| | 15 | 04 | EN-SGP1504-1 |
| | 20 | 04 | EN-SGP2004-1 |
| | 25 | 04 | EN-SGP2504-1 |
| | 30 | 04 | EN-SGP3004-1 |
| | 35 | 04 | EN-SGP3504-1 |
| | 40 | 04 | EN-SGP4004-1 |
| | 45 | 04 | EN-SGP4504-1 |
| | 50 | 04 | EN-SGP5004-1 |
| EdgeFile™ .06 | | | |
| | 15 | 06 | EN-SGP1506-1 |
| | 20 | 06 | EN-SGP2006-1 |
| | 25 | 06 | EN-SGP2506-1 |
| | 30 | 06 | EN-SGP3006-1 |
| | 35 | 06 | EN-SGP3506-1 |
| | 40 | 06 | EN-SGP4006-1 |
| | 45 | 06 | EN-SGP4506-1 |
| | 50 | 06 | EN-SGP5006-1 |
| EdgeOne Fire™ | | | |
| | Small | | EN-EOFSMGP-1 |
| | Primary | | EN-EOFPRIGP-1 |
| | Medium | | EN-EOFMEDGP-1 |
| | Large | | EN-EOFLGGP-1 |
| | Assorted | | EN-EOFGPAST-1 |
| EdgeSequel Sapphire™ | | | |
| 15 | 04 | | EN-SGP1504-1 |
| | 06 | | EN-SGP1506-1 |
| 20 | 04 | | EN-SGP2004-1 |
| | 06 | | EN-SGP2006-1 |
| 25 | 04 | | EN-SGP2504-1 |
| | 06 | | EN-SGP2506-1 |
| 30 | 04 | | EN-SGP3004-1 |
| | 06 | | EN-SGP3006-1 |
| 35 | 04 | | EN-SGP3504-1 |
| | 06 | | EN-SGP3506-1 |
| 40 | 04 | | EN-SGP4004-1 |
| | 06 | | EN-SGP4006-1 |
| 45 | 04 | | EN-SGP4504-1 |
| | 06 | | EN-SGP4506-1 |
| 50 | 04 | | EN-SGP5004-1 |
| | 06 | | EN-SGP5006-1 |

PAPER POINTS - Pack of 100

| Family | Tip Size | Taper Size | SKU |
|---|----------|------------|---------------|
| EdgeTaper™ & EdgeTaper Platinum™ | | | |
| F1 | 20 | | EN-ETPPF1-1 |
| F2 | 25 | | EN-ETPPF2-1 |
| F3 | 30 | | EN-ETPPF3-1 |
| F4 | 40 | | EN-ETPPF4-1 |
| F5 | 50 | | EN-ETPPF5-1 |
| EdgeFile™ .04 | | | |
| | 15 | 04 | EN-SPP1504-1 |
| | 20 | 04 | EN-SPP2004-1 |
| | 25 | 04 | EN-SPP2504-1 |
| | 30 | 04 | EN-SPP3004-1 |
| | 35 | 04 | EN-SPP3504-1 |
| | 40 | 04 | EN-SPP4004-1 |
| | 45 | 04 | EN-SPP4504-1 |
| | 50 | 04 | EN-SPP5004-1 |
| EdgeFile™ .06 | | | |
| | 15 | 06 | EN-SPP1506-1 |
| | 20 | 06 | EN-SPP2006-1 |
| | 25 | 06 | EN-SPP2506-1 |
| | 30 | 06 | EN-SPP3006-1 |
| | 35 | 06 | EN-SPP3506-1 |
| | 40 | 06 | EN-SPP4006-1 |
| | 45 | 06 | EN-SPP4506-1 |
| | 50 | 06 | EN-SPP5006-1 |
| EdgeOne Fire™ | | | |
| | Small | | EN-EOFSMPP-1 |
| | Primary | | EN-EOFPRIIP-1 |
| | Medium | | EN-EOFMEDPP-1 |
| | Large | | EN-EOFLGPP-1 |
| | Assorted | | EN-EOFPAST-1 |
| EdgeSequel Sapphire™ | | | |
| 15 | 04 | | EN-SPP1504-1 |
| | 06 | | EN-SPP1506-1 |
| 20 | 04 | | EN-SPP2004-1 |
| | 06 | | EN-SPP2006-1 |
| 25 | 04 | | EN-SPP2504-1 |
| | 06 | | EN-SPP2506-1 |
| 30 | 04 | | EN-SPP3004-1 |
| | 06 | | EN-SPP3006-1 |
| 35 | 04 | | EN-SPP3504-1 |
| | 06 | | EN-SPP3506-1 |
| 40 | 04 | | EN-SPP4004-1 |
| | 06 | | EN-SPP4006-1 |
| 45 | 04 | | EN-SPP4504-1 |
| | 06 | | EN-SPP4506-1 |
| 50 | 04 | | EN-SPP5004-1 |
| | 06 | | EN-SPP5006-1 |

THE NEW GOLD STANDARD FOR PERFORMANCE, FLEXIBILITY & STRENGTH



EdgeFile X7™

STERILE - PACK OF 6

- Heat-treated FireWire™ NiTi gives amazing flexibility, capable of 90° curves
- 2-8 times more cyclic fatigue resistant than other file systems.*
- Tips 17 - 45 .04 & .06 tapers
- Available in 21mm, 25mm & 29mm
- Max 1mm flute diameter
- Unmatched flexibility
- No shape memory
- Universal compatibility

Gutta Percha - Pack of 60

| Family | Tip Size | Taper Size | SKU |
|-------------------------|----------|------------|----------------|
| EdgeFile™ X7 .04 | | | |
| | 17 | 04 | EN-X71704GP-1 |
| | 20 | 04 | EN-X72004GP-1 |
| | 25 | 04 | EN-X72504GP-1 |
| | 30 | 04 | EN-X73004GP-1 |
| | 35 | 04 | EN-X73504GP-1 |
| | 40 | 04 | EN-X74004GP-1 |
| | 45 | 04 | EN-X74504GP-1 |
| | Assorted | 04 | EN-X7GP04AST-1 |
| EdgeFile™ X7 .06 | | | |
| | 17 | 06 | EN-X71706GP-1 |
| | 20 | 06 | EN-X72006GP-1 |
| | 25 | 06 | EN-X72506GP-1 |
| | 30 | 06 | EN-X73006GP-1 |
| | 35 | 06 | EN-X73506GP-1 |
| | 40 | 06 | EN-X74006GP-1 |
| | 45 | 06 | EN-X74506GP-1 |
| | Assorted | 06 | EN-X7GP06AST-1 |



PAPER POINTS - Pack of 100

| Family | Tip Size | Taper Size | SKU |
|-------------------------|----------|------------|---------------|
| EdgeFile™ X7 .04 | | | |
| | 15 | 04 | EN-SPP1504-1 |
| | 20 | 04 | EN-SPP2004-1 |
| | 25 | 04 | EN-SPP2504-1 |
| | 30 | 04 | EN-SPP3004-1 |
| | 35 | 04 | EN-SPP3504-1 |
| | 40 | 04 | EN-SPP4004-1 |
| | 45 | 04 | EN-SPP4504-1 |
| | 50 | 04 | EN-SPP5004-1 |
| | Assorted | 04 | EN-SPP04AST-1 |
| EdgeFile™ X7 .06 | | | |
| | 15 | 06 | EN-SPP1506-1 |
| | 20 | 06 | EN-SPP2006-1 |
| | 25 | 06 | EN-SPP2506-1 |
| | 30 | 06 | EN-SPP3006-1 |
| | 35 | 06 | EN-SPP3506-1 |
| | 40 | 06 | EN-SPP4006-1 |
| | 45 | 06 | EN-SPP4506-1 |
| | 50 | 06 | EN-SPP5006-1 |
| | Assorted | 06 | EN-SPP06AST-1 |



EdgeSequel Sapphire™

STERILE - PACK OF 4

- Heat-treated FireWire™ NiTi
- 2X more resistant to cyclic fatigue compared to competitor files*
- Triangular cross section with open flute design
- Constant Taper
- ISO Sizes 15-90 .04 tapers
- ISO Sizes 15-50 .06 tapers
- Available in 21mm, 25mm and 31mm



EdgeGlidePath™

STERILE - PACK OF 4

- Heat-treated FireWire™ NiTi
- Replacing as many as 8 instruments to create a faster glidepath
- Replaces #15, 20, 25, 30, 35 hand files in a step-back technique and a #1, 2 & 3 Gates Glidden
- Available in 21mm, 25mm, and 31mm



EdgeTaper Platinum™

STERILE - PACK OF 6

- Heat-treated FireWire™ NiTi gives amazing flexibility, capable of 90° curves
- Familiar technique
- 2X more resistance to cyclic fatigue compared to competitor files*
- SX, S1, S2, F1, F2, F3, F4, F5
- 19mm, 21mm, 25mm, 31mm



EdgeV-Taper HT™

NON STERILE - PACK OF 6

- Heat-treated FireWire™ NiTi
- The Original Minimal preparation file
- Changing Variable Taper
- Tips 20, 22, 25, 30, 35, 40, 45 and 50
- GlidePath 13/.03, 14/.03 and 17/.04
- 21mm, 25mm and 31mm (29mm)



EdgeOne Fire™

STERILE - PACK OF 3

- Heat-treated FireWire™ NiTi
- EdgeOne Fire™ is designed to shape canals in arc reciprocating motion, in a similar way to leading reciprocating files brands
- EdgeOne Fire™ can be used with your existing motor and hand piece reciprocating settings
- Excellent flexibility, capable of 90° curves
- No bounce back to preserve canal anatomy
- Sizes: Small (20), Primary (25), Medium (35) & Large (45)
- Available in 21mm, 25mm and 31mm

